Spider Screw product range

Description	Length	Thread	Head
Spider K1- Ø1,5	6,5 - 8 -10	Self-tapping and Self drilling	Bracket head
Spider K2 -Ø1,9	5 - 6 - 7 - 9 - 11	Self-tapping and Self drilling	Bracket head
Spider C1- Ø1,5	6,5 - 8 - 10	Self-tapping	Bracket head
Spider PIN Ø1,3	8 - 10	Self-tapping	Button head
Smart Screw Ø 1,6	8 - 10	Self-tapping	Button head
Spider SL1-S - Ø1,5	6,5 - 8 -10	Self-tapping and Self drilling	Self-ligating with slot
Spider SL2-S - Ø1,9	6 - 7 - 9 - 11	Self-tapping and Self drilling	Self-ligating with slot
Spider SL1 -Ø1,5	6,5 - 8 -10	Self-tapping and Self drilling	Self-ligating without slot
Spider SL2 -Ø1,9	6 -7 - 8 -9 - 10 -11	Self-tapping and Self drilling	Self-ligating without slot
Spider Regular Plus K2 Ø1,9	7 - 8 - 9 - 10 - 11 - 12 - 13 - 14 -15 -18	Self-tapping and Self drilling	Button head + internal thread
Spider Reg Plus Konic Ø2,0 - 2,3	7 - 8 - 9 - 10 - 11 - 12 - 13 - 14 -15	Self-tapping and Self drilling	Button head + internal thread

Product description

HDC Spider Screw System consists of titanium self-tapping screws with various shapes and sizes for applications in the orthodontic field. It is intended to serve as a fixed anchorage point for the attachment of orthodontic and pre-prosthetics appliances, in order to facilitate the orthodontic movement of teeth. **Spider Screw** has been developed to offer a number of versatile anchorage options capable of immediate loading. Spider Screw has been cleaned, decontaminated, and **sterilized via gamma radiation** and therefore is ready to use. Its associated **accessories** are supplied **NON STERILE** and should be sterilized before use.

Indications for use

The HDC Spider Screw is a threaded titanium dental implant screw. The device is used temporarily with the intention to be removed after orthodontic treatment. **Screws are intended for single use only.**

Insertion sites

MAXILLA	MANDIBLE
Infrazygomatic crest (*)	
Edentulous ridges	Edentulous ridges
Interdicular spaces	Interdicular spaces
Retromolar region	Retromolar region
Alveolar ridges	Sinphisy
Palate	Alveolar ridges

Note: (*) see figure below

Fig.1 – Infrazygomatic crest

Directions for use

HDC Spider Screw comes wrapped in a special double packaging that guarantees its absolute sterile conditions until ready to use. It simplifies as well the transfer of the Spider Screw from its envelope to the placement site. When inserting the Spider Screw, ensure that all the hygiene measures required for invasive surgery are completed.

Step 1

Choose the right Spider Screw based on need by checking the information written on the label on the cardboard box.

Step 2

Once the blister packaging is extracted from its box, remove the Tyvek label and lay the item down on a sterile surface.

Step 3

Collect the Spider Screw from inside its container and place it on a surgical tray.

Step 4**Spider Screw K1 - K2 and SL1-SL2**

Spider Screw K1, K2, SL1 and SL2 can be inserted without any previous drilling depending upon the bone density. It is the responsibility of the professional to determine suitability on a case by case basis before following the procedure. In some cases, especially in the posterior areas of the mandible and palate, pre-drilling is recommended in order to prevent breakage of the Spider Screw. We advise the use of HDC drill FSC-1108 for Spider Screw K1 and SL1; FSC-1309 for Spider Screw K2 and SL2.

Spider Screw C1

Self-tapping Spider Screw C1 requires pre-drilling. We advise the use of HDC drill FSC-1210. In some cases, pre-drilling may be avoided depending on the bone density and based upon the operating physician's recommendation.

Spider Screw PIN

Self-tapping Spider Screw PIN requires pre-drilling. We advise the use of HDC drill FSC-0910. In some cases, it can be inserted without any pre-drilling depending upon the bone density. It is the responsibility of the professional to determine suitability on a case by case basis before following the procedure.

Spider Screw Regular Plus + Konic

Spider Screw Regular Plus can be inserted without any previous drilling depending upon the bone density. It is the responsibility of the professional to determine suitability on a case by case basis before following the procedure. In some cases, especially in the posterior areas of the mandible and palate, pre-drilling is recommended in order to prevent breakage of the Spider Screw. We advise the use of HDC drill FSC-1309.

Step 5

Unscrew the container cap, exposing the Spider Screw in its slot.

Step 6**Spider Screw K1-K2 and C1**

We advise a manual insertion of Spider Screw K1-K2-C1. Pick up the sterile Spider Screw and transfer it to the prepared site for placement using Driver DSP-5052S + DSX-1690S.

Mechanical insertion

It is also possible the use of DPQ-2820 connected to a slow speed handpiece for picking up the sterile Spider Screw K1-K2-C1 and transfer it to the prepared site of placement. The insertion force must not exceed 20 N/cm.

Spider Screw PINManual insertion

Pick up the sterile Spider Screw PIN and transfer it to the prepared site for placement using Driver DSX-1690S + DSP-2352S or DST-1600 + DSP-2352S.

Mechanical insertion

Pick up the sterile Spider Screw and transfer it to the prepared site for placement using Driver DPQ-2322 connected to a slow speed handpiece. The insertion force must not exceed 15 N/cm.

Spider Screw SL1-S and SL2-S

We advise a manual insertion of Spider Screw SL1-SL2 with Slot. Pick up Spider Screw SL1/SL2 S and transfer it to the prepared site for placement using Driver DSX-1690S + DSP-5652S or DST-1600 + DSP-5652S.

Mechanical insertion

It is also possible the use of Driver DPQ-3420 connected to a slow speed handpiece for picking up the sterile Spider Screw and transfer it to the prepared site of placement. The insertion force must not exceed 20 N/cm.

Spider Screw Regular Plus + Konic

We advise a manual insertion of Spider Screw Regular Plus. Pick

up the sterile Spider Screw and transfer it to the prepared site for placement using Driver DSX-1690S + DSP-3852S or DST-1600 + DSP-3852S.

Mechanical insertion

It is also possible the use of Driver DPQ-2820 (Regular Plus) or DPQ-3825 (Regular Plus Konic) connected to a slow speed handpiece for picking up the sterile Spider Screw and transfer it to the prepared site of placement. The insertion force must not exceed 20 N/cm.

Step 7**Spider Screw K1-K2-C1**

Use Handle Driver DSQ-2824 or Long Handle Driver DSX-1690S along with Driver DSX-2852S for final positioning tuning.

Spider Screw PIN


Use Handle Driver DSQ-2324 for final positioning tuning.

Spider Screw SL1-S and SL2-S

Use Handle Driver DSQ-3424 for SL1-S and SL2-S for final positioning tuning.

Step 8**Spider Screw SL1-SL2**

After the bone screw insertion, its head appears like in fig. 1. It is possible to use orthodontic wires .018"x.022 or .018"x.025 (fig. 2). To do this, insert the device into the Spider Screw head's slot and turn it clockwise using DXL-2820 Hand Drive device, until the wire is fixed (fig. 3). Maximum rotation angle 90°.

**Step 9****Load Resistance**

Spider Screw has been developed to offer a number of versatile anchorage options capable of immediate loading. Immediate loading is possible because the Spider Screw is a non-osteointegrable implant and consequently force can be applied immediately after placement. When loading, please be careful to apply a perpendicular load (90°) with reference to the long axis of the Spider Screw. The applied force can range from 50g to 300g depending on screw choice, bone quality, and the desired orthodontic movement.

Contraindications

- Insufficient bone quantity and/or poor bone quality in the receiving site;
- Poor oral hygiene;
- Heavy smoking, tobacco and alcohol abuse;
- Systemic blood disorders;
- Uncontrolled diabetes.

Warnings

- Placement of surgical screws requires specific knowledge of anatomy and techniques and this procedure must be carried out by qualified and trained people. Improper patient selection and/or incorrect technique can cause placement failure and/or loss of supporting bone.
- An implanted device or used screw should never be reused. Any screws which have been contaminated with blood or body fluids should be discarded.

Precautions

- Effective and complete screening of screw application candidate must be performed. Visual inspection as well as panoramic and periapical radiographs are essential to determine anatomical landmarks and bone adequacy. Lateral or other types of X-ray examinations are recommended.
- Detailed instructions, limitations, and possible adverse effects of the procedure should be given to the patient.
- Spider Screw application surgical procedures have some risks which include the injury of delicate anatomical structures both of the superior jaw-bone and of mandibular one, if existing conditions

are not carefully considered.

- **HDC Spider Screw** has been designed to achieve a flexible anchorage with immediate loading and of limited duration. Consequently, the efficiency of this system should be not depending from osteointegration process. A trained and qualified personnel should be aware that HDC Spider Screw is not designed as for osteointegrable anchorage (deferred loading).

Adverse effects

After Spider Screw application, untimely anchorage loss may occur. Potential causes include but are not limited to:

- Bone poor quantity and/or quality, Osteoporosis, Osteolysis, Osteomyelitis, inhibited revascularization, or infection can cause loosening, bending, cracking or fracture of the device or premature loss of fixation with the bone leading to non-union.
 - Infections.
 - Poor oral hygiene or patient's cooperation and/or genetic diseases (diabetes).
 - Migration, bending, fracture or loosening of the implant.
 - Metal sensitivity or allergic reaction to a foreign body.
 - Pain, discomfort, or abnormal sensation due to the presence of the device.
 - Increased fibrous tissue response around the fracture site and/or the implant.
 - Necrosis of bone.
 - Inadequate healing.
 - Localized swelling, oedema, and tissue reaction.
- Apart from these adverse effects, there are always possible complications of any surgical procedure such as, but not limited to, infection, nerve damage, and pain which may not be related to the implant.

Shipping and handling

Every **HDC Spider Screw®** has been cleaned, decontaminated, and sterilized via gamma radiation and, therefore, is ready to use. Sterile packages containing the screw should be opened upon a sterile surface. See "Directions for use".

Caution

This medical device has not to be reused.

Sterility guaranteed unless package is opened or damaged.

Sales of this medical device are subjected to provisions of directive 2007/47/CEE of 05/09/20074 and of D.L. 37 of 25/01/2010.

FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

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HDC Srl